


Centrifuges



Model	Ultra 8 Variable	Ultra 8 Select	Ultra 8 Digital	E8 Fixed	E8 Variable	Portafuge
Model #	ULC-08AV-1501	ULC-08AS-1501	ULC-08AD-1501	E8C-08AF-1501	E8C-08AV-1501	E8C-08AF-150P
Type	Test Tube	Test Tube	Test Tube	Test Tube	Test Tube	Test Tube
Features	Variable Speed	4 Speed Selectable	Programmable	Fixed Speed	Variable Speed	Fixed Speed
Capacity	8 place X 15ml	8 place X 15ml	8 place X 15ml	8 place X 15ml	8 place X 15ml	8 place X 15ml
Tube Size Min/Max	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 17 X 130mm
Max g	1400	1400	1400	1400	1400	1500
Max Speed (rpm)	3300	3300	3300	3300	3300	+/-3300
Uses	Blood, Urine, Semen, Fecal	Blood, Urine, Semen, Fecal	Blood, Urine, Semen, Fecal	Blood Only	Blood, Urine, Semen, Fecal	Blood Only



Model	C3 Select	C5	C5 - 24 Place	CxR
Model #	C3C-08SS-1501	C5C-08SV-15T1	C5C-BKSV-15T3	CXC-06SF-1501
Type	Test Tube Swing Out	Test Tube Swing Out	Test Tube	Test Tube Swing Out
Features	Selectable Pre-set Speeds to 3300	5000 RPM Programmable	5000 RPM Programmable	3400 RPM Short Profile
Capacity	8 place X 15ml *4 place X 50ml	8 place X 15ml *4 place X 50ml	24 X 2.5-10ml *4 place X 50ml tube inserts	6 place X 15ml
Tube Size Min/Max	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 16 X 130mm	10.25 X 65mm / 17 X 130mm	10.25 X 65mm / 16 X 130mm
Max g	1500	4000	4000	1800
Max Speed (rpm)	3000	5000	5000	3400
Uses	Blood, Urine, Semen, Fecal	Blood (STAT), Urine, Semen, Fecal, Hematocrits, PRP/PPP, Stem-Cell, Irap	Blood (STAT), Urine, Semen, Fecal, Hematocrits, PRP/PPP, Stem-Cell	Blood Only



*Optional Features

P 770.270.1394
F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30046

LWScientific.com

Centrifuges



Model	V24 Combo	V24 Combo Tach	M24 Hematocrit	Zipocrit	ZipSpin
Model #	CMC-24AV-1501 (choose rotor)	CMC-24AV-15T1 (choose rotor)	HEC-24HF-7501	ZPC-04HF-750P	ZPC-06AF-020P
Type	Combination	Combination	Microhematocrit	Microhematocrit	Microcentrifuge
Features	Test Tube & Microhematocrit	Test Tube & Microhematocrit	Heavy Duty 24 Place	Portable angled hematocrit	Portable 6 & 12 place rotors
Capacity	24 X 15ml or 24 X 75mm	24 X 15ml or 24 X 75mm	24 X 75mm	4 X 75mm	(6 & 12) X 2ml
Tube Size Min/Max	10.25 X 65mm /17X130mm or 75mm capillary	10.25X65mm/17X130mm or 75mm capillary	75mm capillary (only)	75mm capillary (only)	0.5ml / 2.0ml
Max g	2450 - 14900	2450 - 14900	14900	4400	2800
Max Speed (rpm)	4000 - 12000	4000 - 12000	4000 - 12000	12000	11000
Users	Blood, Urine, Semen, Fecal, Hematocrits, DNA Separation	Blood, Urine, Semen, Fecal, Hematocrits, DNA Separation	Hematocrits Only	Hematocrits Only	Blood, Urine, DNA Separation



P 770.270.1394
F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30046

LWScientific.com