

VARIABLE LOUPES

LPM-P35V-4407

Zoom magnification **2.5x to 3.5x**
Variable working distance 440mm-250mm
Variable depth of field 130mm-50mm
Variable field of view 100mm-40mm



Weight 2.5oz



3.5X PREMIUM LOUPES



Weight 2.5oz

SHORT: Close-up, bench work

LPM-P35F-4007 (3.5x)

-Working distance 365 +/- 15mm
-Depth of field 30mm
-Field of view 44mm

LONG: Stand-up surgeries

LPM-P35F-4607 (3.5x)

-Working distance 460 +/- 20mm
-Depth of field 40mm
-Field of view 45mm

2.5X PREMIUM LOUPES

SHORT: Close-up, bench work

LPM-P25F-3007 (2.5x)

-Working distance 340 +/- 25mm
-Depth of field 50mm
-Field of view 66mm

LONG: Stand-up surgeries

LPM-P25F-4607 (2.5x)

-Working distance 460 +/- 30mm
-Depth of field 80mm
-Field of view 90mm



Weight 2.6oz



LPM-C25F-3807 (2.5x)

-Working distance 380 +/- 30mm
-Depth of field 60mm
-Field of view 80mm

LPM-C35F-3807 (3.5x)

-Working distance 380 +/- 20mm
-Depth of field 40mm
-Field of view 50mm

Works on most prescription glasses

CLIP-ON LOUPES



Weight 2.3 oz

